

DISCOVER THE

# CATALYST *ADVANTAGE*



CATALYST 5



CATALYST 5Ti



CATALYST 5X



CATALYST 4



# CATALYST


Since 2006, Catalyst has been the best performing and most durable ultra-lightweight folding chair on the market. Catalyst continues to provide one-of-a-kind features, enhancements and exclusive innovations to the complex rehab industry.

Catalyst is designed to make your team more efficient. Every Catalyst model is constructed using the same high-quality components and innovative design.

← → Folded Width - 4"




## BETTER *OPTIONS & ACCESSORIES*



### **SUPER LOW SWING AWAY**

Lower the client's seat by 1.5 inches, making this the perfect option for foot propulsion.



### **HANGER RELEASE CHOICES**

Three options to meet the client's or caregivers needs. Choose from a classic hanger, 4-way inline or 4-way out.



### **X-HINGE DUAL CROSS BRACE**

The dual interlocking cross brace provides the rigidity needed to maximize the client's efficiency.



### **MAXX PERFORMANCE**

The superior rigidity, performance and responsiveness gives the client more out of every push.




### **PRO ELR**

Enhanced adjustability in the pivot point gives the client a personalized fit that can improve clinical outcomes.



### **RESIDUAL LIMB SUPPORT**

The ease of adjustment and swing in swing out design allows the client to achieve the optimal angle and position with ease and without tools.



### **PERFORMANCE CASTER HOUSING**

Easily adjust the client's caster housing to achieve their optimal setup with an additional 1-1.5" of seat slope.



### **TRANSIT WHEELS**

Offer the client the ability to navigate narrow spaces and doorways with rear frame integrated Transit Wheels.

## **AXIOM SEATING**

Experience Axiom's vast offering of backrests, cushions and custom seating that can be installed on Catalyst. Get simplified ordering through the ACS App.



# CATALYST *ADVANTAGE*

K0005 Comparables*	Ki Mobility Catalyst 5	Ki Mobility Catalyst 5 VX	Sunrise Medical Quickie 2	Sunrise Medical QXi	TiLite Aero X	Motion Composites Helio A7	Motion Composites Helio A6
Standard Weight Capacity	300 lbs	300 lbs	300 lbs	300 lbs	265 lbs	265 lbs	265 lbs
Standard Weight †	23 lbs	24 lbs	25.4 lbs	28.8 lbs	25.2 lbs	Not Published	Not Published
Transport Weight §	12.8 lbs	13.9 lbs	13.6 lbs	17.6 lbs	18 lbs	13.6 lbs	13 lbs
Base Retail (USD)	\$2,120	\$1,895	\$2,157	\$1,902	\$2,390	\$2,295	\$2,095
Cross Brace Design	Dual w/ Power X-Hinge	Dual w/ Power X-Hinge	Dual	3-Arm	3-Arm	Single	Single
Lifetime Warranty on Footrest	•	•					
Frame Material	7000 Series Aluminum	7000 Series Aluminum	Aluminum	Aluminum	Aluminum	7000 Series Aluminum	6061 Series Aluminum
Titanium Frame Upgrade	•				•		
WC-19 Ansi-Resna Transit Approved	•	•	•	•	•	•	•
Make It Right Program	•	•	•	•			
Free Growth Program	•		•				
Heavy Duty Option ( 350 lbs )	•	•	•		•	•	•

\* All information published on manufacturers website as of 2/18/2021

† 16x16 standard configuration, including rear wheels and footrests

§ 16x16 standard configuration, without rear wheels and footrests

