

Jazzy[®] EVO 613 SERIES

THE MOST ADVANCED JAZZY[®] EVER

You will love the all-new Medicare-reimbursable Jazzy[®] EVO 613 series! You can choose between the EVO 613 or the EVO 613 Li, which uses a lithium iron phosphate battery for longer life, more convenient charging and fewer service calls.



Scan QR code
to learn more.



25%
FASTER
CHARGING TIME

55%
LIGHTER
& SMALLER BATTERY

300%
LONGER
BATTERY LIFE CYCLE

13"
WHEELS

3"
GROUND
CLEARANCE



Get the lithium
advantage with the
Jazzy[®] EVO 613 Li.

FDA Class II
Medical Device¹

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

800.800.8586

pridemobility.com/Jazzy-Evo-613-Li



EVO 613 SERIES

AVAILABLE MODELS

Jazzy® EVO 613 / Jazzy® EVO 613 Li
Group 2, Standard (K0822/K0823)



FDA Class II
Medical Device¹

Black (matte)

White (matte)

Red (matte)

Robin's Egg Blue
(matte)

Iceberg Blue
(matte)

Sugar Plum
(matte)

**Lithium Iron
Phosphate Battery**
(Available on Jazzy® EVO 613 Li)

- Faster charging time
- Sustained power performance
- Lighter weight
- More useable energy
- Longer life



FEATURES

- Maximum speed up to 5.1 mph
- 13" mid-wheel drive tires
- Narrow base width of only 22"
- Group 2, Standard (K0822/K0823)
- Choice of six contemporary removable, matte-finish color-through shrouds
- Depth adjustable, limited-recline, memory foam, high-back with single seat post allows the user to swivel
- Pride's patented Active-Trac® (ATX) Suspension
- Mid-Wheel 6® Technology

ACCESSORIES

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag kit
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- 2.1A XLR USB charger

OPTIONS

- Elevated leg rests
- Swing-away leg rests
- Residual limb support
- Heel loops
- Unoccupied transit kit
- Depth-adjustable solid seat pan, limited recline, high-back with swivel

SPECIFICATIONS

WEIGHT CAPACITY 300 lbs.

MAXIMUM SPEED² Up to 5.1 mph

DIMENSIONS³ Width: 22" Length: 41" with footrest

TOTAL WEIGHT EVO 613: 189 lbs. EVO 613 Li: 160.5 lbs.

TURNING RADIUS 22"

RANGE PER CHARGE³ EVO 613: Up to 14 miles @ 300 lbs.
EVO 613: Up to 17.1 miles @ 200 lbs.
EVO 613 Li: Up to 17.6 miles @ 300 lbs.
EVO 613 Li: Up to 21 miles @ 200 lbs.

WARRANTY

FRAME: 5-year limited
DRIVETRAIN: 2-year limited
ELECTRONICS: 2-year limited

SLA BATTERY: 13-month limited
LITHIUM BATTERY: 3-year limited

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

AUTHORIZED PROVIDER



pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]