

# Durable medical equipment **product guide**












Innovative products and support



# Elevate care with medical equipment that goes the distance

Find lots of new products inside including folding options for bath safety items and additional arm and seat options for wheelchairs.

## What's inside:

	Mobility rack.....	3
	Bath safety.....	4
	Wheelchairs.....	14
	Transport chairs.....	22
	Rollators.....	25
	Walkers.....	33
	Knee walkers.....	38
	Homecare beds & accessories.....	39
	Canes.....	42

# Mobility rack

Free with \$3,500 minimum purchase order of Durable Medical Equipment products

## Recommended planogram:

### Walking Aids

Item No.	Description	Pkg.	Page
G07850M	Flip tray for Guardian walkers	3/cs	37
MDS86615K	Walker basket	2/cs	37
MDS86410KDBW	Upright folding walker, 5" wheels	1/ea	35
MDS4PTCANETIP	4-point cane tips	4/cs	43
MDS86425W	Quad cane tips	6/cs	43
MDS86426W4	Cane tips	4/cs	43
MDS86222W	Quad cane, small base, Black	2/cs	42
MDS86420	Offset cane, Black	6/cs	42
MDS86228	Quad cane, large base, Black	2/cs	42
MDS86430W	Folding T-handle cane	6/cs	42
4UCANE	Folding hybrid cane, offset, Black	1 ea	43
4UCANEB	Folding hybrid cane, offset, Blue	1 ea	43
MDS86850ESKD	Basic steel rollator, Burgundy	1 ea	27
MDS86615SKI	Ski gliders for walkers	6 pr/cs	37
MDS86615T	Walker tips	1 pr	37
MDSWBALL	Walker ball glides	6 pr/cs	37
MDS86615BAG	2-button walker bag	4/cs	37
MDSRUBSKI	Walker glides	6/cs	37
MDS86615L	10" leg replacement	1 pr	37
MDS86615W5	Footpiece, 5" wheels	1 pr	37
MDS86615W5T	5" long extension wheels	1 pr	37
MDS808200KDR	Lightweight steel transport chair	1 ea	22

### Aids for Daily Living

Item No.	Description	Pkg.	Page
MDS86031RCE	31" reacher, Blue	6/cs	40
MDS6800BA	Bed assist bar, Black	3/cs	40
MDS6800SWIVRH	SwivAssist bar	1 ea	40

### Bath Safety

Item No.	Description	Pkg.	Page
G2-101KRX1	Bath chair with back	1 ea	7
G2-201KRX1	Bath chair without back	1 ea	8
G3-100KBRX1	Transfer bench with back	1 ea	9

MDSMOBRACK



### 3-in-1 commodes

- Use as a bedside commode or over toilet
- Folding commode reduces storage space by 35%

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	HCPCS	Retail Pkg.	Pkg.
MDS89664	Grey, fixed, standard seating	350 lbs. (158.8 kg)	15"-21"	15"	19"	24"	E0163	N	4/cs
G1-200X1	White, folding, standard seating	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	E0163	N	1 ea
G1-200X4	White, folding, standard seating	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	E0163	N	4/cs
G1-200RX1	White, folding standard seating	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	E0163	Y	1 ea
G1-200LX1	White, folding, elongated seating	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	E0163	N	1 ea
G1-200LX4	White, folding, elongated seating	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	E0163	N	4/cs
G1-200MX1	Dark Grey folding standard seating with Microban*	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	-	N	1 ea
G1-200LMX1	Dark Grey folding elongated seating with Microban*	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	-	N	1 ea

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

### Drop arm commodes

- Armrests conveniently swing away for easy patient transfer
- Padded options available for increased comfort

Drop arm holds up to 350 lbs.

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	HCPCS	Retail Pkg.	Pkg.
G1-301DX1	White, drop arm, elongated	350 lbs. (158.8 kg)	19"-23"	16.5"	20"	23"	E0165	N	1 ea
G1-301DX2	White, drop arm, elongated	350 lbs. (158.8 kg)	19"-23"	16.5"	20"	23"	E0165	N	2/cs
G1-301DPX1	White, padded, drop arm, elongated	350 lbs. (158.8 kg)	18"-22"	17"	20"	23"	E0165	N	1 ea
G1-303DWX1	White drop arm elongated, 4 locking casters	300 lbs. (136 kg)	17.5"-21.5"	16.5"	17.5"	23.5"	-	N	1 ea
G98204	Chrome padded seat and back with drop arm	350 lbs. (158.8 kg)	20"-25"	18"-20"	18.5"	22"	E0165	N	1 ea



G1-200X1



G1-200LX1



G1-301DX1



G1-303DWX1

Easy to fold  
and budget  
friendly

## Economy commodes

- Basic commode without seat lid or plastic armrests
- Available with standard or elongated seat

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	HCPCS	Retail Pkg.	Pkg.
G1-202EX4	White, folding standard seating	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	E0163	N	4/cs
G1-202ELX4	White, folding elongated seating	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	E0163	N	4/cs

## Bariatric commodes

- Fixed and folding commodes have sturdy steel frames
- Perfect for high weight capacity needs

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	HCPCS	Retail Pkg.	Pkg.
G1-201LBX1	White, bariatric steel folding commode	650 lbs. (294.8 kg)	16"-22"	16.5"	24"	31"	E0168	N	1 ea
G1-400BX2	White, bariatric steel fixed bench seat commode	650 lbs. (294.8 kg)	17"-20.5"	18.5"	24"	29.5"	E0168	N	2/cs
G1-400BDX1	White, bariatric steel drop arm commode	1,000 lbs. (453.6 kg)	17.5"-21.5"	18.5"	25"	31"	E0168	N	1 ea



G1-202ELX4



G1-201LBX1



G1-400BX2



G1-400BDX1

## Aluminum commodes

- Rust resistant aluminum
- Multi-functional bathroom usage

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	HPCPS	Retail Pkg.	Pkg.
G1-501LWX1	Aluminum elongated commode, with four locking casters	300 lbs. (136.1 kg)	21.5"	16.5"	18"	22"	-	N	1 ea
G1-502PWX1	Aluminum padded shower chair commode, with four locking casters	300 lbs. (136.1 kg)	20.5"	16"	18.5"	21"	E0163	N	1 ea
G1-503WCPX1	Aluminum shower commode wheelchair	275 lbs. (124.7 kg)	20"-22"	15.5"	19"	25"	-	N	1 ea
G1-504WCPX1	Aluminum shower commode transport chair	275 lbs. (124.7 kg)	20"-22"	15.5"	19"	23"	-	N	1 ea
G1-505LX1	Aluminum 3-in-1 folding commode	350 lbs. (158.8 kg)	15.5"-21.5"	16.5"	17"	22"	-	N	1 ea
G1-506DWX1	Aluminum drop arm commode, with four locking casters	300 lbs. (136.1 kg)	18.5"-23"	16.5"	17.5"	22"	-	N	1 ea
G1-507BWX1	Aluminum bariatric commode, with four locking casters	1,000 lbs. (453.6 kg)	21"	21"	25.5"	28.5"	-	N	1 ea



G1-502PWX1



G1-503WCPX1



G1-507BWX1

## Buckets

- Replacement commode buckets
- Available with or without lid and handle

Item No.	Description	HPCPS	Retail Pkg.	Pkg.
MDS80306B	Grey bucket with lid and metal handle	E0167	N	6/cs
MDS80306EB	Grey bucket only	E0167	N	18/cs
G1-004X18	Grey bucket lid only	-	N	18/cs
G1-005X6	Grey bucket with lid and plastic handle	E0167	N	6/cs



MDS80306B

## Commode liners

- Disposable commode liners make clean up easy

Item No.	Description	HPCPS	Retail Pkg.	Pkg.
MDS89664LINER	Commode liner with absorbent pad, 12 liners per roll	-	Y	6 rolls/cs
MDS89664LINEH	Commode liner with absorbent pad, 12 liners per roll	-	Y	1 roll
G1-002X1	Commode liners, no absorbent pad, 50 liners per each	-	N	1 ea



MDS89664LINER

## Shower chairs with back

- Rust resistant aluminum chairs are perfect for patients who have difficulty standing in the shower
- Backrest included for added stability, security and comfort

Item No.	Description	Wt. capacity	Seat height	Seat depth	HCPCS	Retail Pkg.	Pkg.
G2-100KH1	Grey shower chair with arms	350 lbs. (158.8 kg)	16"-21"	13.5"	-	N	1 ea
MDS89745RAH	White shower chair with arms	350 lbs. (158.8 kg)	16"-21"	13.5"	E2045	Y	1 ea
MDS89745ELMBH	Tan shower chair with arms, Microban*	350 lbs. (158.8 kg)	15"-21"	15"	-	Y	1 ea
MDS89745KDMBH	Light blue shower chair, Microban*	300 lbs. (136.1 kg)	13.75"-20.75"	12"	-	Y	1 ea
G2-101KX1	Grey tool-free shower chair	400 lbs. (181.4 kg)	13.75"-20.5"	12"	E0240	N	1 ea
G2-101KR1	Grey retail shower chair, tool-free	400 lbs. (181.4 kg)	13.75"-20.5"	12"	E0240	Y	1 ea
G2-102KX1	Grey knockdown shower chair	400 lbs. (181.4 kg)	15"-19"	12"	E0240	N	1 ea
G2-102KX4	Grey knockdown shower chair	400 lbs. (181.4 kg)	15"-19"	12"	E0240	N	4/cs
G2-103RX1	Grey folding tool-free shower chair	300 lbs. (136.1 kg)	14.5"-20.5"	15"	-	Y	1 ea
G2-104AX2	Grey deluxe aluminum shower chair with back	400 lbs. (181.4 kg)	13.5"-20.5"	12"	-	N	2/cs
G2-105RX1	White premium plastic shower chair with back	350 lbs. (158.8 kg)	16"-20.5"	15"	-	Y	1 ea
MDSMOMCHAIRG	Dark Grey Momentum shower chair	250 lbs. (113.4 kg)	16"-20.5"	15"	-	Y	2/cs
MDSMOMCHAIRGH	Dark Grey Momentum shower chair	250 lbs. (113.4 kg)	16"-20.5"	15"	-	Y	1 ea
G2-111KM1	White shower chair, Microban*	400 lbs. (181.4 kg)	16.5"-20.5"	14"	-	N	1 ea
G2-112KM1	Blue shower chair, Microban*	400 lbs. (181.4 kg)	16.5"-20.5"	14"	-	N	1 ea

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



G2-100KH1



G2-102KX1



G2-105RX1

## Shower chairs without back

- Rust resistant aluminum shower chairs for patients who have difficulty standing in the shower
- No backrest for easy maneuverability in tight spaces

Item No.	Description	Wt. capacity	Seat height	Seat depth	HCPCS	Retail Pkg.	Pkg.
G2-200KHX1	Grey shower chair with arms	350 lbs. (158.8 kg)	16" -21"	13.5"	-	N	1 ea
MDS89740RWAH	White shower chair with arms	350 lbs. (158.8 kg)	16"-21"	13.5"	-	Y	1 ea
G2-201KRX1	Grey retail shower chair, tool-free	400 lbs. (181.4 kg)	13.75"-20.5"	12"	E0245	Y	1 ea
G2-202KX1	Grey shower chair, knockdown	400 lbs. (181.4 kg)	15"-19"	12"	E0245	N	1 ea
G2-203FRX1	Grey folding shower chair	250 lbs. (113.4 kg)	14.5"-16"	12"	E0245	Y	1 ea
G2-204AX2	Grey assembled shower chair	400 lbs. (181.4 kg)	13.5"-20.5"	12"	-	N	2/cs
G2-300KRX1	Grey bath stool	300 lbs. (136.1 kg)	13.75"-20.5"	12.5"	E0245	Y	1 ea
G2-300KRPX1	Grey padded bath stool with swivel	300 lbs. (136.1 kg)	16"-22"	13.75"	E0245	Y	1 ea
G2-400RX1	Grey tub portable bath bench	250 lbs. (113.4 kg)	22"	12.5"	E0245	Y	1 ea



G2-200KHX1



G2-201KRX1



G3-300KRPX1

Great for small spaces



### Bariatric shower chairs

- Heavy duty reinforced frames
- For high weight capacity needs

Item No.	Description	Wt. capacity	Seat height	Seat depth	HCPCS	Retail Pkg.	Pkg.
G2-100BAX1	Grey assembled shower chair with back	650 lbs. (294.8 kg)	16.5"-20.5"	17"	E0245	N	1 ea
G2-102BX1	Grey shower chair with back	550 lbs. (249.5 kg)	15"-20"	12"	E0245	N	1 ea
G2-202BX1	Grey shower chair without back	550 lbs. (249.5 kg)	15"-20"	12"	E0245	N	1 ea

### Transfer benches

- Allows patients to enter the bath tub from a seated position
- Padded and bariatric options available to meet various needs

Item No.	Description	Wt. capacity	Seat height	Seat depth	HCPCS	Retail Pkg.	Pkg.
G3-100KBX1	Grey transfer bench	400 lbs. (181.4 kg)	18"-22.5"	16.25"	E0247	N	1 ea
G3-100KBX2	Grey transfer bench	400 lbs. (181.4 kg)	18"-22.5"	16.25"	E0247	N	2/cs
G3-100KBRX1	Grey transfer bench	400 lbs. (181.4 kg)	18"-22.5"	16.25"	E0247	Y	1 ea
G3-101KBX1	Grey transfer bench commode	400 lbs. (181.4 kg)	18"-22"	16"	E0247	N	1 ea
G3-101KBPX1	Grey padded transfer bench commode	400 lbs. (181.4 kg)	18"-22"	16"	E0247	N	1 ea
G3-200KBX1	Grey bariatric transfer bench	500 lbs. (226.8 kg)	19.25"-23.5"	16.5"	E0248	N	1 ea
MDS86952XW	White bariatric transfer bench	550 lbs. (249.5 kg)	16"-20"	16"	E0248	N	1 ea
G3-001RX1	White transfer bench shower curtain	N/A	N/A	N/A	-	Y	1 ea

550 lbs. weight cap at a fraction of the price



### Economy raised toilet seat

- Makes using the toilet easier for those with trouble bending
- Simple design made of polypropylene resin is easy to clean
- Slip-resistant pads and lip inside toilet offer additional stability

Item No.	Description	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
G4-400RX1	White raised toilet seat with side locks, 4"	250 lbs. (113.4 kg)	-	Y	1 ea
MDS80314	White elevated locking toilet seat, without arms, 5"	400 lbs. (181.4 kg)	E0244	Y	3/cs
MDS80314H	White elevated locking toilet seat, without arms, 5"	400 lbs. (181.4 kg)	E0244	Y	1 ea
MDS80315EL	White elongated elevated locking toilet seat with arms, 5"	350 lbs. (158.8 kg)	-	Y	3/cs
MDS80316	White elevated locking toilet seat, with arms, 5"	350 lbs. (158.8 kg)	E0244	Y	3/cs
MDS80316H	White elevated locking toilet seat, with arms, 5"	350 lbs. (158.8 kg)	E0244	Y	1 ea
MDS80316MB	Light Blue elevated locking toilet seat, with arms, with Microban*, 5"	350 lbs. (158.8 kg)	-	Y	3/cs
G30250	White toilet seat riser, without arms, without lock, 5"	250 lbs. (113.4 kg)	E0244	Y	3/cs
G30250H	White toilet seat riser, without arms, without lock, 5"	250 lbs. (113.4 kg)	E0244	Y	1 ea
G30260	White toilet seat riser, without arms, with lock, 5"	350 lbs. (158.8 kg)	E0244	Y	3/cs
G30270A	White toilet seat riser, with arms, with lock, 5"	350 lbs. (158.8 kg)	E0244	Y	3/cs
G30270AH	White toilet seat riser, with arms, with lock, 5"	350 lbs. (158.8 kg)	E0244	Y	1 ea
G4-111MX1	White elevated locking toilet seat, with arms, with Microban*, 5"	400 lbs. (181.4 kg)	-	N	1 ea
G4-112MX1	Blue elevated locking toilet seat, with arms, with Microban*, 5"	400 lbs. (181.4 kg)	-	N	1 ea

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



G4-400RX1



MDS80315EL



MDS80316MB



G30260



G30270A

## White enamel grab bars

- Enamel coated steel grab bars mount to the wall
- Can be used horizontally or vertically and include mounting hardware

Item No.	Dimensions	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
MDS86012	12" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86012H	12" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G87520M	12" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86016	16" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86016H	16" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G87530M	16" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86018	18" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86018H	18" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G87540M	18" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86024	24" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
G87550M	24" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86032	32" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
G87560M	32" wide	300 lbs. (136.1 kg)	-	Y	3/cs
G87570M	36" wide	300 lbs. (136.1 kg)	-	Y	3/cs



## Knurled chrome grab bars

- Chrome plated steel grab bars mount to the wall
- Can be used horizontally or vertically and include mounting hardware

Item No.	Dimensions	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
MDS86012CHR	12" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86012CHRH	12" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G81020KEM	12" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86016CHR	16" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86016CHRH	16" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G81030KEM	16" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86018CHR	18" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86018CHRH	18" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G81040KEM	18" wide	300 lbs. (136.1 kg)	-	Y	3/cs
MDS86024CHR	24" wide	250 lbs. (113.4 kg)	E0241	Y	3/cs
MDS86024CHRH	24" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
MDS86032CHR	32" wide	250 lbs. (113.4 kg)	-	Y	3/cs
MDS86032CHRH	32" wide	250 lbs. (113.4 kg)	E0241	Y	1 ea
G81050KEM	24" wide	300 lbs. (136.1 kg)	E0241	Y	3/cs
G81060KEM	32" wide	300 lbs. (136.1 kg)	-	Y	3/cs
G81070KEM	36" wide	300 lbs. (136.1 kg)	-	Y	3/cs



### Suction grab bars

- Attach to the wall without screws
- Grab bars stay in place with suction cups

Item No.	Description	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
MDS86012SCLI	White with locking indicator, 12"L	250 lbs. (113.4 kg)	-	Y	3/cs
MDS86016SCLI	White with locking indicator, 16"L	250 lbs. (113.4 kg)	-	Y	3/cs
MDS86019SCLI	White, pivot grab bar, 19"L	250 lbs. (113.4 kg)	-	Y	3/cs
MDS86012SCLIH	White with locking indicator, 12"L	250 lbs. (113.4 kg)	-	Y	1 ea
MDS86016SCLIH	White with locking indicator, 16"L	250 lbs. (113.4 kg)	-	Y	1 ea
MDS86019SCLIH	White, pivot grab bar, 19"L	250 lbs. (113.4 kg)	-	Y	1 ea



MDS86012SCLI

### Tub grab bars

- Step-through clamp design for better hand positioning
- Tool-free installation with lined pads that won't mar tub surf

Item No.	Description	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
MDS86321W	White, locks to side of tub	250 lbs. (113.4 kg)	E0246	Y	2/cs
MDS86321WH	White, locks to side of tub	250 lbs. (113.4 kg)	E0246	Y	1 ea
G98006	White, bi-level grab bar	250 lbs. (113.4 kg)	E0241	Y	3/cs



MDS86321W

### Cast protectors

- Help keep casts dry when bathing
- Available in adult and child sizes

Item No.	Description	HCPCS	Retail Pkg.	Pkg.
CUR200ALL	Clear, adult leg	-	Y	6/cs
CUR200ALLH	Clear, adult leg	-	Y	1 ea
CUR200AAA	Clear, adult arm	-	Y	6/cs
CUR200AAAH	Clear, adult arm	-	Y	1 ea
CUR100KLL	Clear, child leg	-	Y	6/cs
CUR100KLLH	Clear, child leg	-	Y	1 ea
CUR100KAA	Clear, child arm	-	Y	6/cs
CUR100KAAH	Clear, child arm	-	Y	1 ea



CUR200ALL



CUR200AAA

### Shower heads

- Easy to install white handheld shower heads with diverter valve
- Two styles

Item No.	Description	HCPCS	Retail Pkg.	Pkg.
G7-100RX1	Deluxe hand held shower 3 massage functions	-	Y	1 ea
MDSHSHOWER	White hand held with 6.5 ft. hose	-	Y	6/cs
MDSHSHOWERH	White hand held with 6.5 ft. hose	-	Y	1 ea



G7-100RX1

### Momentum shower caddies

- Basket helps organize shampoo, conditioner and soap
- Supports up to 5 lbs

Item No.	Description	HCPCS	Retail Pkg.	Pkg.
MDSMOMCADDYG	Wall mounted shower caddy	-	Y	2/cs
MDSMOMCADDYGH	Wall mounted shower caddy	-	Y	1 ea



MDSMOMCADDYG

### Bath mats

- Slip-resistant heavy weight mat with suction cups on bottom for secure grip
- Antimicrobial treatment helps protect bath mat\* from growth of stains and odors caused by bacteria, mold and mildew

Item No.	Description	HCPCS	Retail Pkg.	Pkg.
MDSBATHMAT	White bath mat with Microban*	-	Y	4/cs
MDSBATHMATH	White bath mat with Microban*	-	Y	1 ea
MDSMOMMATG	Grey Momentum bath mat with Microban*	-	Y	2/cs
MDSMOMMATGH	Grey Momentum bath mat with Microban*	-	Y	1 ea

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



MDSMOMMATGH

### Bath furniture

- Teak look without the constant maintenance
- Frames are aluminum seat and shelf are composite wood

Item No.	Description	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
G2-500X1	Bathroom bench with faux wood	200 lbs. (90.7 kg)	-	N	1 ea
G2-501X1	Bathroom corner bench with faux wood	200 lbs. (90.7 kg)	-	N	1 ea



G2-500X1



G2-501X1

# Wheelchairs

## K1 Guardian wheelchairs

- Durable hammertone frame; available in nylon or vinyl upholstery
- Hemi-height adjustable
- Accommodates 2" seat extender

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length swing back	Swing-away	K1166N22S	Black	300 lbs. (136.1 kg)	16" x 16"	K001 or K002	N	1 ea
		K1186N22S	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea
		K1206N22S	Black	300 lbs. (136.1 kg)	20" x 16"	K001 or K002	N	1 ea
	Elevating	K1166N22E	Black	300 lbs. (136.1 kg)	16" x 16"	K001 or K002	N	1 ea
		K1186N22E	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea
		K1206N22E	Black	300 lbs. (136.1 kg)	20" x 16"	K001 or K002	N	1 ea
Full-length removable	Swing-away	K1186N13S	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea
	Elevating	K1186N13E	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea

### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length swing back	Swing-away	K1166V22S	Black	300 lbs. (136.1 kg)	16" x 16"	K001 or K002	N	1 ea
		K1186V22S	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea
		K1206V22S	Black	300 lbs. (136.1 kg)	20" x 16"	K001 or K002	N	1 ea
	Elevating	K1166V22E	Black	300 lbs. (136.1 kg)	16" x 16"	K001 or K002	N	1 ea
		K1186V22E	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea
		K1206V22E	Black	300 lbs. (136.1 kg)	20" x 16"	K001 or K002	N	1 ea
Full-length removable	Swing-away	K1186V11S	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea
	Elevating	K1186V11E	Black	300 lbs. (136.1 kg)	18" x 16"	K001 or K002	N	1 ea

Swing-away  
arm rest

Accommodates  
2" seat depth  
extenders



K1166V22S

## K2 Guardian wheelchairs

- Durable hammertone frame
- Ships in hemi height; includes seat belt, anti tippers and brake extenders

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Full-length removable and height adjustable	Elevating	K2166N42E1	Black	300 lbs. (136.1 kg)	16" x 16"	K0001 or K0002+E0971, E0973, E0978, E0961, K0195	N	1 ea
		K2186N42E1	Black	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002+E0971, E0973, E0978, E0961, K0196	N	1 ea
		K2206N42E1	Black	300 lbs. (136.1 kg)	20" x 16"	K0002 + E0973 + K0195	N	1 ea



K2166N42E1

## K4 lightweight wheelchairs

- Durable, lightweight TIG-welded frame folds down for easy transport

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, swing back	Swing-away	MDS806500N	Black	300 lbs. (136.1 kg)	16" x 16"	K0004	N	1 ea
		MDS806500	Black	300 lbs. (136.1 kg)	18" x 16"	K0004	N	1 ea
		MDS806560	Black	350 lbs. (158.8 kg)	20" x 16"	K0004+E2201	N	1 ea
		MDS806570	Black	350 lbs. (158.8 kg)	22" x 16"	K0004+E2201	N	1 ea
	Elevating	MDS806550N	Black	300 lbs. (136.1 kg)	16" x 16"	K0004+K0195	N	1 ea
		MDS806550	Black	300 lbs. (136.1 kg)	18" x 16"	K0004+E2201	N	1 ea
		MDS806565	Black	350 lbs. (158.8 kg)	20" x 16"	K0004+E2201+K0195	N	1 ea
		MDS806575	Black	350 lbs. (158.8 kg)	22" x 16"	K0004+E2201+K0195	N	1 ea
Full-length removable	Swing-away	MDS806500FLA	Black	300 lbs. (136.1 kg)	18" x 16"	-	N	1 ea
	Elevating	MDS806550FLA	Black	300 lbs. (136.1 kg)	18" x 16"	-	N	1 ea



MDS806560

### K4 Basic wheelchairs

- Durable grey frame
- Hemi height adjustable

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, swing back	Swing-away	MDS806500NE	Silver	300 lbs. (136.1 kg)	16" x 16"	K0004	N	1 ea
		MDS806500E	Silver	300 lbs. (136.1 kg)	18" x 16"	K0004	N	1 ea
		MDS806500NE18	Silver	300 lbs. (136.1 kg)	16" x 18"	K0004	N	1 ea
		MDS806560E	Silver	300 lbs. (136.1 kg)	20" x 18"	K0006 or K0003+E2201	N	1 ea
	Elevating	MDS806550NE	Silver	300 lbs. (136.1 kg)	16" x 16"	K0004+K0195	N	1 ea
		MDS806550E	Silver	300 lbs. (136.1 kg)	18" x 16"	K0004+K0195	N	1 ea
		MDS806500E18	Silver	300 lbs. (136.1 kg)	18" x 18"	-		1 ea
		MDS806565E	Silver	300 lbs. (136.1 kg)	20" x 18"	K0006 or K0003+E2201+K0195	N	1 ea

### K4 Plus wheelchairs

- Durable hammertone frame
- Includes anti-tippers

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, swing-back and height adjustable	Swing-away	MDS806500PLUS	Black	300 lbs. (136.1 kg)	18" x 16"	-	N	1 ea
		MDS806560PLUS	Black	300 lbs. (136.1 kg)	20" x 18"	K0004+E2201+E0971	N	1 ea
	Elevating	MDS806550PLUS	Black	300 lbs. (136.1 kg)	18" x 16"	-	N	1 ea



MDS806500E



MDS806550PLUS



### K3 Guardian

- Durable hammertone frame and easy-to-clean nylon upholstery
- Swing-back arms are also removable for easy transfers and storage
- Hemi-height adjustable

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HPCPS	Retail Pkg.	Pkg.
Desk-length swing-back & removable	Swing-away	K3166N24S	Black	300 lbs. (136.1 kg)	16" x 16"	K0003	N	1 ea
		K3186N24S	Black	300 lbs. (136.1 kg)	18" x 16"	K0003	N	1 ea
		K3206N24S	Black	300 lbs. (136.1 kg)	20" x 16"	K0003	N	1 ea
	Elevating	K3166N24E	Black	300 lbs. (136.1 kg)	16" x 16"	K0003+K0095	N	1 ea
		K3186N24E	Black	300 lbs. (136.1 kg)	18" x 16"	K0003+K0095	N	1 ea
		K3206N24E	Black	300 lbs. (136.1 kg)	20" x 16"	K0003+K0095	N	1 ea
Desk-length, height-adjustable, swing-back & removable	Swing-away	K3166N34S	Black	300 lbs. (136.1 kg)	16" x 16"	K0003+E0951+E0971+E0973	N	1 ea
		K3186N34S	Black	300 lbs. (136.1 kg)	18" x 16"	K0003+E0951+E0971+E0973	N	1 ea
		K3206N34S	Black	300 lbs. (136.1 kg)	20" x 16"	K0003+E0951+E0971+E0973	N	1 ea
	Elevating	K3166N34E	Black	300 lbs. (136.1 kg)	16" x 16"	K0003+E0951+E0971+E0974+E0990	N	1 ea
		K3186N34E	Black	300 lbs. (136.1 kg)	18" x 16"	K0003+E0951+E0971+E0974+E0990	N	1 ea
		K3206N34E	Black	300 lbs. (136.1 kg)	20" x 16"	K0003+E0951+E0971+E0974+E0990	N	1 ea
Full-length, swing back & removable	Swing-away	K3166N14S	Black	300 lbs. (136.1 kg)	16" x 16"	K0003	N	1 ea
		K3186N14S	Black	300 lbs. (136.1 kg)	18" x 16"	K0003, K0003	N	1 ea
		K3206N14S	Black	300 lbs. (136.1 kg)	20" x 16"	K0003, K0007+K0195+E0990	N	1 ea
	Elevating	K3166N14E	Black	300 lbs. (136.1 kg)	16" x 16"	K0003, K0003+E0990+K0195	N	1 ea
		K3186N14E	Black	300 lbs. (136.1 kg)	18" x 16"	K0003, K0007+K0195+E0990	N	1 ea
		K3206N14E	Black	300 lbs. (136.1 kg)	20" x 16"	K0003, K0003	N	1 ea
Full-length, height-adjustable, swing back & removable	Swing-away	K3166N44S	Black	300 lbs. (136.1 kg)	16" x 16"	K0003	N	1 ea
		K3186N44S	Black	300 lbs. (136.1 kg)	18" x 16"	K0003	N	1 ea
		K3206N44S	Black	300 lbs. (136.1 kg)	20" x 16"	K0003, K0003	N	1 ea
	Elevating	K3166N44E	Black	300 lbs. (136.1 kg)	16" x 16"	K0003, K0003 + E0990 or K0195	N	1 ea
		K3186N44E	Black	300 lbs. (136.1 kg)	18" x 16"	K0003, K0003 + E0990 or K0195	N	1 ea
		K3206N44E	Black	300 lbs. (136.1 kg)	20" x 16"	K0003, K0003 + E0990 or K0195	N	1 ea



K3186N34S



K3186N14S



K3186N14E

## Guardian recliners

- Infinite position reclining up to 180°
- Durable hammertone frame

## Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Elevating	KR168N21E	Black	300 lbs. (136.1 kg)	16" x 18"	K0006+E0990+E1226	N	1 ea
		KR188N21E	Black	300 lbs. (136.1 kg)	18" x 18"	K0006+E0990+E1226	N	1 ea
		KR208N21E	Black	300 lbs. (136.1 kg)	20" x 18"	K0006+E0990+E1226	N	1 ea
Full-length removable	Elevating	KR188N11E	Black	300 lbs. (136.1 kg)	18" x 18"	K0006+E0990+E1226	N	1 ea

## Reclining wheelchairs

- Infinite position reclining from 90°–140°
- Durable carbon steel frame
- Hemi-height adjustable

## Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Elevating	MDS808350	Chrome	300 lbs. (136.1 kg)	16" x 17"	K0003+E0966+E1225+K0195+E0971	N	1 ea
		MDS808450	Chrome	300 lbs. (136.1 kg)	18" x 17"	K0003+E0966+E1225+K0195+E0972	N	1 ea
		MDS808550	Chrome	300 lbs. (136.1 kg)	20" x 17"	K0003+E0966+E1225+K0195+E0973+E2201	N	1 ea
		MDS808650	Chrome	300 lbs. (136.1 kg)	22" x 17"	K0007 + E1226 + E2201 + E0971 + E0966 + K0195 or E0990	N	1 ea

## Hybrid wheelchairs

- Combination wheelchair and transport chair
- Quick release large wheel

## Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Swing-away	MDS806250NH2	Blue	300 lbs. (136.1 kg)	16" x 16"	K0004 or E1038	N	1 ea
		MDS806250H2	Blue	300 lbs. (136.1 kg)	18" x 16"	K0004 or E1038	N	1 ea
	Elevating	MDS806300H2	Blue	300 lbs. (136.1 kg)	18" x 16"	K0004 or E1038+K0195 or E0990	N	1 ea

KR168N21E



MDS808550



MDS806250H2



## Excel wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Swing-away	MDS806250N	Chrome	300 lbs. (136.1 kg)	16" x 16"	K0001 or K0002	N	1 ea
		MDS806250D	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002 + E0990 or K0195	N	1 ea
		MDS806250NVY	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002 + K0195 or E0990	N	1 ea
		MDS806400	Chrome	300 lbs. (136.1 kg)	20" x 16"	K0001 + E2201	N	1 ea
	Elevating	MDS806300N	Chrome	300 lbs. (136.1 kg)	16" x 16"	K0001 or K0002 + K0195	N	1 ea
		MDS806300D	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002	N	1 ea
		MDS806300NVY	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001	N	1 ea
		MDS806450	Chrome	300 lbs. (136.1 kg)	20" x 16"	K0001 + E2201 + K0195	N	1 ea
Full-length permanent	Swing-away	MDS806150D	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002	N	1 ea
	Elevating	MDS806200N	Chrome	300 lbs. (136.1 kg)	16" x 16"	K0002 + K0195 or E0990	N	1 ea
		MDS806100D	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002	N	1 ea
Full-length removable	Swing-away	MDS806400FLA	Chrome	300 lbs. (136.1 kg)	20" x 16"	-	N	1 ea
	Elevating	MDS806300DFLA	Chrome	300 lbs. (136.1 kg)	18" x 16"	K0001 or K0002	N	1 ea

## Excel extra-wide wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Swing-away	MDS806700	Chrome	500 lbs. (226.8 kg)	20" x 18"	K0007	N	1 ea
		MDS806800	Chrome	500 lbs. (226.8 kg)	22" x 18"	K0007	N	1 ea
		MDS806900	Chrome	500 lbs. (226.8 kg)	24" x 18"	K0007+E2202	N	1 ea
	Elevating	MDS806750	Chrome	500 lbs. (226.8 kg)	20" x 18"	K0007+K0195	N	1 ea
		MDS806850	Chrome	500 lbs. (226.8 kg)	22" x 18"	K0007+K0195	N	1 ea
		MDS806950	Chrome	500 lbs. (226.8 kg)	24" x 18"	K0007+E2202+K0195	N	1 ea
Full-length removable	Swing-away	MDS806700FLA	Chrome	500 lbs. (226.8 kg)	20" x 18"	K0007	N	1 ea
		MDS806900FLA	Chrome	500 lbs. (226.8 kg)	24" x 18"	K0007+E2202	N	1 ea
		MDS806750FLA	Chrome	500 lbs. (226.8 kg)	20" x 18"	K0007+K0195	N	1 ea
	Elevating	MDS806850FLA	Chrome	500 lbs. (226.8 kg)	22" x 18"	K0007+K0195	N	1 ea
		MDS806950FLA	Chrome	500 lbs. (226.8 kg)	24" x 18"	K0007+E2202+K0195	N	1 ea



MDS806300D



MDS806900

## Shuttle bariatric wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Elevating	MDS809650	Black	850 lbs. (385.5 kg)	24" x 20"	K0007+K0108+E2203+K0195	N	1 ea
		MDS809750	Black	850 lbs. (385.5 kg)	26" x 20"	K0007+K0108+E2203+K0195	N	1 ea
		MDS809850	Black	850 lbs. (385.5 kg)	28" x 20"	K0007+K0108+E2203+K0195	N	1 ea
		MDS809950	Black	850 lbs. (385.5 kg)	30" x 20"	K0007+K0108+E2203+K0195	N	1 ea

## Pediatric wheelchair

- Adjustable telescoping handles.
- Hemi-height adjustable

### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, swing-back	Elevating	MDS806140PD	Yellow	250 lbs. (113.4 kg)	14" x 12"	E1238+E0971+K0195+E0990	N	1 ea

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, swing-back	Swing-away	KPD2N22S	Yellow	250 lbs. (113.4 kg)	12" x 12"	-	N	1 ea
		KPD4N22S	Yellow	250 lbs. (113.4 kg)	14" x 12"	-	N	1 ea
		KPD2N22S1	Blue	250 lbs. (113.4 kg)	12" x 12"	-	N	1 ea
		KPD4N22S1	Yellow	250 lbs. (113.4 kg)	14" x 12"	-	N	1 ea



MDS809750



KPD2N22S



KPD4N22S1



MDS806140PD

## Wheelchair accessories

Item No.	Description	Pkg.
MDS85190	Holder, oxygen tank, IV pole, wheelchair	1 ea
WCA6100	Holder, oxygen tank, for wheelchair ADL	1 ea
MDS85189	Device, anti-tip, wheelchair, rear	1 pr
STMSM22	Anti-rollback device	1 ea
MDS85189H	Bag, leg rest or carry, all for wheelchairs	1 pr
MDS85181U	Universal oxygen tank holder for wheelchairs	1 ea
MDS85183	Wheelchair IV pole, telescoping	1 ea
MDS85189WK4	Anti-tip rear wheel for K4 and recliner	1 pr
WCABELT	Safety belt for wheelchair	1 ea
MDS85157PR	Black handle extension wheel lock 5/8 x 8	1 pr
MDS85197A	Adjustable anti-theft device for wheelchair	1 ea
MDSK41RMDRIVE	One arm drive attachment for Excel K	1 ea
WCA806991	Nut, tinnerman for wheelchair, leg rigging	1 pr
WCA6200	Cup holder, ADL	1 ea
MDS85196	Anti-fold device	1 ea
MDS85195	Kit, super hemi, with standard wheel	1 ea
MDS851965SH	Anti-fold device for shuttle wheelchair	1 ea
MDS85197SH	Anti-theft, overhead for shuttle wheelchair	1 ea
STMSM22A	Anti-rollback device with alarm	1 ea
MDS85180	5-in-1 IV/O2/anti-theft device	1 ea
STMSM22W	Anti-rollback device	1 ea
MDS85189SH	Anti-tip, bariatric, for shuttle wheelchair	1 pr
MDS85189K	Anti-tip Kidz chair	1 pr
WCAA0060	K1/K2 nylon seat extender, 18"	1 ea
WCAA0061	K1/K2 nylon seat extender, 20"	1 ea
WCAA0070	K1 & K2 anti-tippers	1 pr
WCAA0069	K3 anti-tipper	1 pr
WCAA0066	K3 16" seat extender	1 ea
WCAA0067	K3 18" seat extender	1 ea
WCAA0068	K3 20" seat extender	1 ea



MDS85183

# Transport chairs

## Basic steel transport chairs with 8" wheels

- Durable steel frame with black upholstery
- Includes detachable footrests and a seat belt for safety
- Back folds down for convenient transporting



MDS808200B

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	HCPCS	Retail Pkg.	Pkg.
MDS808200B	Black	250 lbs. (113.4 kg)	Swing-away	Full-length permanent	19"	E0138	N	1 ea
MDS808200BN	Black	250 lbs. (113.4 kg)	Swing-away	Full-length permanent	17"	E0138	N	1 ea
MDS808200WAT	Red	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	-	N	1 ea
MDS808200	Chrome	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	E0138	N	1 ea
MDS808200KDR	Burgundy	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	-	Y	1 ea



MDS808200ABE

## Aluminum transport chairs with 8" wheels

- Powder coated aluminum frame with black upholstery
- Includes seat belt for safety and padded arms

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	HCPCS	Retail Pkg.	Pkg.
MDS808200ABE	Blue	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	E1038	N	1 ea
MDS808200ARE	Red	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	E1038	N	1 ea

### Aluminum transport chairs with 12" wheels

- Powder coated aluminum frame with black upholstery
- Includes seat belt for safety and padded arms
- Larger back tire for maneuvering through obstacles
- Hand brakes on handles



MDS808210ARE

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	HCPCS	Retail Pkg.	Pkg.
MDS808210ABE	Blue	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	E1038	N	1 ea
MDS808210ARE	Red	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	E1038	N	1 ea

### Excel® Translator combination rollator/transport chairs

- Durable steel frame with black upholstery
- Includes easy-to-adjust, push-button, swing-away footrests
- Footrests lock to sides of rollator when not in use
- Height-adjustable push grips and comfortable hand brake



MDS808200TR

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	HCPCS	Retail Pkg.	Pkg.
MDS808200TRR	Red	250 lbs. (113.4 kg)	Swing-away	Desk length permanent	19"	-	N	1 ea
MDS808200TR	Blue	250 lbs. (113.4 kg)	Swing-away	Desk length permanent	19"	E0143+ E0156	N	1 ea

### Ultralight transport chairs

- Lightweight aluminum frame with black upholstery
- Carrying handles and compact design allow for easy and convenient storage and transport



MDS808200F3B

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	HCPCS	Retail Pkg.	Pkg.
MDS808200F3R	Red	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	-	Y	1 ea
MDS808200F3B	Blue	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	-	Y	1 ea
MDS808200F3BK	Black	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	-	N	1 ea
MDS808200F3S	Silver	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	-	N	1 ea

### Bariatric transport chairs

- Durable frame with black upholstery, seat belt, handbrakes and wheel lock for better control



MDS80822B

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	HCPCS	Retail Pkg.	Pkg.
TCB228N21S	Bariatric Steel, Blue	500 lbs. (226.8 kg)	Swing-away	Desk length, swing back	22"	-	N	1 ea
MDS80822B	Bariatric steel, Red	500 lbs. (226.8 kg)	Swing-away	Desk-length	22"	E1039	N	1 ea
MDS808200BAR	Bariatric, lightweight, aluminum, Red	400 lbs. (181.4 kg)	Elevating	Elevating	22"	E1039	N	1 ea



# Rollators

## Basic steel rollators

- Durable powder-coated steel frame
- Storage bag under the seat is included

Item No.	Description	Wt. capacity	Seat height	Wheel size	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86850ESKD	Burgundy	350 lbs. (158.8 kg)	20.5"	6"	31"-35"	17"	E0143+E0156	Y	1 ea
MDS86850ES	Burgundy	350 lbs. (158.8 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86850EBS	Blue	300 lbs. (136.1 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86830EBSL	Black	350 lbs. (158.8 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86850ES8	Burgundy	350 lbs. (158.8 kg)	20.5"	8"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86850EBS8	Blue	300 lbs. (136.1 kg)	20.5"	8"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86830EBSL8	Black	300 lbs. (136.1 kg)	20.5"	8"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86840EBS	Light Blue	300 lbs. (136.1 kg)	20.5"	6"	31"-36"	19.75"	E0143+E0156	N	1 ea
MDS86840EBS8	Light Blue	300 lbs. (136.1 kg)	20.5"	8"	31"-36"	19.75"	E0143+E0156	N	1 ea



MDS86840EBS



MDS86850E8

### Basic aluminum rollators

- Durable aluminum frame
- Storage bag under the seat is included



MDS86850EB

Item No.	Description	Wt. capacity	Seat height	Wheel size	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86850EB	Blue	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86850E	Burgundy	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	Y	1 ea
MDS86850EP	Purple	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86850EG	Green	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86830EBL	Black	300 lbs. (136.1 kg)	20.5"	6"	31"-35"	17.5"	E0143+E0156	N	1 ea
MDS86850EJR	Burgundy	250 lbs. (113.4 kg)	20.5"	6"	29"-33"	17.5"	E0143+E0156	N	1 ea
MDS86830EBL8	Black	350 lbs. (158.8 kg)	20.5"	8"	31"-35"	17.5"	E0143+E0156	N	1 ea

Aluminum options weigh only 11 lbs!



MDS86825SLBLK

### Superlight rollators with 6" wheels

- Lightweight aluminum frame
- Weighs only 11 lbs. making it easy to transport or store when not in use

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86825SLR	Burgundy	250 lbs. (113.4 kg)	20.5"	29.5"-36"	17"	E0143+E0156	Y	1 ea
MDS86825SLBLK	Blue	250 lbs. (113.4 kg)	20.5"	29.5"-36"	17"	E0143+E0156	Y	1 ea
MDS86825SLB	Purple	250 lbs. (113.4 kg)	20.5"	29.5"-36"	17"	E0143+E0156	Y	1 ea
MDS86825SLRS	Red	350 lbs. (158.8 kg)	18"-23"	29.5"-36"	17"	-	Y	1 ea

\*MDS86825SLRS made of steel

### Bariatric rollators with 8" wheels

- Heavy-duty bariatric rollators
- Wider distance between handles for users needing extra room

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86800XWS	Steel, Burgundy	500 lbs. (226.8 kg)	20.5"	31.5"–37.5"	23.5"	E0149+E0156	N	1 ea
MDS86800XW	Aluminum, Burgundy	400 lbs. (181.4 kg)	20.5"	31.5"–37.5"	23.5"	E0149+E0156	N	1 ea

### Guardian rollators with 8" wheels

- Aluminum frame
- Padded seat and curved padded backrest provide enhanced support and comfort

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
G07887B	Blue	300 lbs. (136.1 kg)	20"	31"–36.5"	17.5"	E0143+E0156	N	1 ea
G07887R	Red	300 lbs. (136.1 kg)	20"	31"–36.5"	17.5"	E0143+E0156	N	1 ea

### Steel rollators with basket

- Easily removable wired basket
- Flip up backrest for easy storage

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86860ERS8	Red	300 lbs. (136.1 kg)	21.25"	30.5"–36"	18.5"	–	N	1 ea
MDS86860EBS8	Blue	300 lbs. (136.1 kg)	21.25"	30.5"–36"	18.5"	–	N	1 ea



### Momentum rollators with 6" wheels

- Lightweight steel frame with large padded seat
- Seat flips up to enable user to walk inside
- All-terrain wheels make it easier to get around
- Product weighs 18.8 lbs.

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86870R	Red	250 lbs. (113.4 kg)	20.5"	34.25"–41.5"	18"	–	N	1 ea
MDS86870W	White	250 lbs. (113.4 kg)	20.5"	34.25"–41.5"	18"	–	N	1 ea
MDS86870G	Grey	250 lbs. (113.4 kg)	20.5"	34.25"–41.5"	18"	–		1 ea

### Empower rollators with 8" wheels

- Aluminum frame with memory foam seat
- Padded backrest flips forward and backward enabling the user to sit comfortably on either side
- Seat flips up to enable user to walk inside

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86845BM	Blue	300 lbs. (136.1 kg)	20.5"	34.25"–41.5"	18"	–	Y	1 ea
MDS86845BLKM	Black	300 lbs. (136.1 kg)	20.5"	34.25"–41.5"	18"	–	Y	1 ea
MDS86845SM	Silver	300 lbs. (136.1 kg)	20.5"	34.25"–41.5"	18"	–	Y	1 ea

### Three wheeled rollators

- Easy to maneuver with a tray, basket and storage bag

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86840TRIB	Blue, 8" wheels	275 lbs. (124.7 kg)	20.5"	31.75"–35.5"	14.5"	–	N	1 ea
MDS86840TRIB	Red, 8" wheels	275 lbs. (124.7 kg)	20.5"	31.75"–35.5"	14.5"	–	N	1 ea



MDS86870G

Seat flips up to allow user to walk inside the walker



MDS86845BM

Includes reflectors on the frame



MDS86840TRIB

## European style rollators

Item No.	Description	Wt. capacity	Seat height	Height range	Color	HCPCS	Retail Pkg.	Pkg.
MDS86890EUB	Simplicity 2 rollator	300 lbs. (136.1 kg)	19"	32"-36"	Blue	-	N	1 ea
MDS86890EUS	Simplicity 2 rollator	300 lbs. (136.1 kg)	19"	32"-36"	Silver	-	N	1 ea
MDS86870ANB	Lightweight knock down	300 lbs. (136.1 kg)	21"	33.8"-38"	Navy Blue	-	N	1 ea
MDS86870SHW	Aluminum with suspension	300 lbs. (136.1 kg)	22"	34.9"-39"	White	-	N	1 ea
MDS86880R	Aluminum euro style	300 lbs. (136.1 kg)	21"	30"-35"	Red	-	N	1 ea
MDSUPROLLW	Upright rollator	300 lbs. (136.1 kg)	21.5"	40.75"-45.5"	White	-	N	1 ea



MDS86890EUB



MDS86870ANB



MDS86870SHW



MDS86880R



MDSUPROLLW

## byACRE carbon fiber rollators

- World's lightest rollators at 10.6 lbs.
- Carbon fiber frame is 5x stronger than steel
- Ergonomic handles have integrated brake lines and are height adjustable in 1" increments
- All byACRE rollators include organizer bag

World's  
lightest  
rollator

### Compact

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	HCPCS	Retail Pkg.	Pkg.
BYA100SMBK	Black	285 lbs. (129.3 kg)	22"	31"-36"	16"	-	N	1 ea
BYA100SMR	Red	285 lbs. (129.3 kg)	22"	31"-36"	16"	-	N	1 ea
BYA100SMW	White	285 lbs. (129.3 kg)	22"	31"-36"	16"	-	N	1 ea

### Regular

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	HCPCS	Retail Pkg.	Pkg.
BYA100MDBK	Black	285 lbs. (129.3 kg)	24"	31"-36"	16.5"	-	N	1 ea
BYA100MDR	Red	285 lbs. (129.3 kg)	24"	31"-36"	16.5"	-	N	1 ea
BYA100MDW	White	285 lbs. (129.3 kg)	24"	31"-36"	16.5"	-	N	1 ea

### Comfort

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	HCPCS	Retail Pkg.	Pkg.
BYA100LGBK	Black	285 lbs. (129.3 kg)	24"	31"-36"	18.5"	-	N	1 ea
BYA100LGR	Red	285 lbs. (129.3 kg)	24"	31"-36"	18.5"	-	N	1 ea
BYA100LGW	White	285 lbs. (129.3 kg)	24"	31"-36"	18.5"	-	N	1 ea



BYA100SMBK



BYA100MDW



BYA100LGR

Folds for easy  
portability



### byACRE all terrain carbon fiber rollator

- Air filled tires
- Weighs 14.78 lbs.
- Carbon fiber frame is 5x stronger than steel
- Ergonomic handles have integrated brake lines and are height adjustable in 1' increments
- Includes organizer bag



BYA1500G

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	HCPCS	Retail Pkg.	Pkg.
BYA1500G	Defender Green	330 lbs (149.7 kg)	24"	32.3"-37.4	18.5"	-	N	1 ea

### byACRE rollator accessories

Item No.	Description	Retail Pkg.	Pkg.
BYAGRBG	Grocery bag, Grey	N	1 ea
BYAWKNBG	Weekend bag, Grey	N	1 ea
BYATVLBG	Travel bag, Grey	N	1 ea
BYACNHR	Cane holder, right	N	1 ea
BYACNHL	Cane holder, left	N	1 ea
BYAOBGR	Organizer bag, small/medium, Grey	N	1 ea
BYAOBGL	Organizer bag, large, Grey	N	1 ea
BYASMBR	Backrest, small/medium	N	1 ea
BYALBR	Backrest, large	N	1 ea



BYAGRBG



BYAWKNBG



BYATVLBG



BYAOBGR



BYASMBR

## Rollator accessories

Item No.	Description	Retail Pkg.	Pkg.
MDS86825BSK	Basket for MDS86825 rollator	N	1 ea
MDS86810BSK	Basket for MDS86810 series rollators	N	1 ea
MDS86800XWBSK	Basket for MDS86800XW rollator	N	1 ea
MDS86850BAG	Black bag for MDS86850 rollators	N	1 ea
MDS86835SBAG	Zebra stripe bag for under seat	N	1 ea
MDS86835WBAG	Black bag for under seat	N	1 ea
MDS86845CH	Cup holder for empower rollators	N	1 ea
MDS86835SPG	Rollator front bag, Grey	N	1 ea
MDS86835SPT	Front rollator bag, Teal	N	1 ea
MDS86835SPR	Front rollator bag, Burgundy	N	1 ea
MDESCUPCANEHW	Combo, cup and cane holder	N	1 ea
MDS86835SM	Mirror	N	1 ea





# Walkers

## Hemi-walkers

- Lightweight design allows the hemi sidestepper to be maneuvered easily—even with one arm
- Bi-level hand grip offers extra assistance for users

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
G07770X1	Adult hemi walker	250 lbs. (113.4 kg)	PVC	32"–36"	N/A	E0135	N	1 ea
G07770XW	Bariatric adult hemi walker	500 lbs. (226.8 kg)	PVC	32"–36"	N/A	E0148	N	1 ea
G07771	Youth hemi walker	250 lbs. (113.4 kg)	PVC	28.5"–32"	N/A	E0135	N	2/cs

## Two-button folding walkers\*

- Durable walker with two-button folding system
- Push-button, height-adjustable legs

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS864104	Aluminum, no wheels	300 lbs. (136.1 kg)	PVC	32"–39"	18"	E0135	N	4/cs
MDS864104B	Steel, no wheels	400 lbs. (181.4 kg)	PVC	32"–39"	17.5"	E0135	N	4/cs
MDS86410W4	Aluminum, 3" wheels	300 lbs. (136.1 kg)	PVC	32"–38"	17.5"	E0143	N	4/cs
MDS86410W54	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	32"–38"	17"	E0143	N	4/cs
MDS86410W54B	Steel, 5" wheels	400 lbs. (181.4 kg)	PVC	32.5"–39.5"	17"	E0143 or E0149	N	4/cs



G07770X1



MDS864104



MDS86410W54B

\*All MDS walkers come unassembled

### Youth two-button folding walkers

- Youth size durable walker with two-button folding system
- Comfortable grip with nonmarking tips



MDS86410JW54

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86410J4	Aluminum, no wheels	300 lbs. (136.1 kg)	PVC	25"-32"	17.5"	E0135	N	4/cs
MDS86410J4B	Steel, no wheels	400 lbs. (181.4 kg)	PVC	25"-32"	17.5"	-	N	4/cs
MDS86410JW54	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	25"-32"	17.5"	E0143	N	4/cs
MDS86410JW54B	Steel, 5" wheels	400 lbs. (181.4 kg)	PVC	25"-32"	17.5"	E0149	N	4/cs

### Bariatric folding walkers

- Extra-wide frames and higher weight capacities provide a comfortable fit for larger patients
- Includes two-button folding system
- All bariatric folding walkers come assembled



MDS86410XWW

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86410XW	Aluminum, no wheels	500 lbs. (226.8 kg)	PVC	31"-41"	20.5"	E0148	N	1 ea
MDS86410XWB	Steel, no wheels	650 lbs. (294.8 kg)	PVC	31"-41"	20"	-	N	1 ea
MDS86410XWW	Aluminum, 5" wheels	500 lbs. (226.8 kg)	PVC	33.5"-43.25"	20.5"	E0149	N	1 ea
MDS86410XWWB	Steel, 5" wheels	600 lbs. (272.2 kg)	PVC	31"-41"	20"	E0149	N	1 ea

## Folding paddle walkers

- Ideal design for patients with manual dexterity limitations
- The paddle mechanism can be flipped-up or pushed-down to fold the walker

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86410P	Aluminum, no wheels	300 lbs. (136.1 kg)	PVC	32"-39"	16.5"	E0135	N	4/cs
MDS86410P5	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	32.5"-39.5"	16.5"	E0143	N	4/cs
MDS86410P5J	Junior, 5" wheels	300 lbs. (136.1 kg)	PVC	25.5"-32.5"	16.5"	E0143	N	4/cs
MDS86410KDBW	Aluminum, 5" wheels	350 lbs. (158.8 kg)	PVC	32.5"-39.5"	16.75"	E0149	Y	1 ea

## Upright folding walkers

- Upright folding walkers help users with posture compared to standard walkers
- Comfortable grip with mark-resistant tips

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDSUPWALK	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	46"-53"	20"	-	N	1 ea
MDSUPWALKR	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	46"-53"	20"	-	N	1 ea
MDSUPWALKJR	Aluminum, youth 5" wheels	300 lbs. (136.1 kg)	Cushion	35.5"-42.5"	18"	-	N	1 ea



MDS86410P



MDS86410KDBW



MDSUPWALK

\*All MDS walkers come unassembled

MDS86410TRIG



### Folding trigger walker

- Trigger mechanism can be pulled up to fold the walker
- Wheeled glidebrakes to help user brake easily
- Walker folds and comes with its own bag for easy storage when travelling
- Two storage pouches included

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
MDS86410TRIG	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	29.5"–37"	17"	–	N	1 ea

### Guardian walkers

- All Guardian walkers come assembled
- Easy-fold sides for easy movement through narrow spaces
- Comfortable grip and mark-resistant tips

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	HCPCS	Retail Pkg.	Pkg.
G30757W	Aluminum, 5" wheels, adult	300 lbs. (136.1 kg)	PVC	32.5"–39.5"	16.5"	E0143	N	4/cs
G30758W	Aluminum, 5" wheels, youth	300 lbs. (136.1 kg)	PVC	27.25"–34.25"	16.5"	E0143	N	4/cs
G30754B	Steel, no wheels, adult	650 lbs. (294.8 kg)	PVC	32"–39"	22"	E0148	N	2/cs
G30755P	Aluminum, no wheels, adult	300 lbs. (136.1 kg)	PVC	32.5"–39.5"	16.5"	E0135	N	4/cs
G30756P	Aluminum, no wheels, youth	300 lbs. (136.1 kg)	PVC	27"–34"	16.5"	E0135	N	4/cs
G07783M	Aluminum, 5" wheels, pediatric	175 lbs (79.4 kg)	PVC	24"–36"	9.5"	E0135	N	1 ea
G07749	Aluminum, no wheels, youth	200 lbs. (90.7 kg)	PVC	24.5"–28.5"	12.75"	E0135	N	1 ea
G07750	Aluminum, no wheels	175 lbs (79.4 kg)	PVC	20.5"–24.5"	10.75"	E0135	N	1 ea
G07767	Aluminum, no wheels, adult	400 lbs. (181.4 kg)	PVC	32.5"–39.25"	19"	E0148	N	1 ea



G30754B



G30755P



G07749

## Walker accessories

Item No.	Description	Retail Pkg.	Pkg.
G07729M	7/8" footpiece, 3" wheels and glide caps	Y	1 pr
G07731-8	1" footpiece, 3" wheels and glide caps	Y	1 pr
G07731G	Footpiece, 3" wheels	Y	3 pr/cs
G07722-8	Guardian, 5" front walker wheels	Y	1 pr
MDS86615W	3" casters for walkers	Y	2 pr/cs
MDS86615W5	Footpiece, 5" wheels and glide caps	Y	1 pr
MDS86615W5W	Footpiece, 5" wheels and glide caps	Y	6 pr/cs
G07724	Pediatric extension legs, 5" wheels	Y	2 pr/cs
G07725	Walker footpiece extension kit, 5" wheels	Y	1 kit
G07825-8	Swivel footpiece, 3" wheels	Y	1 pr
G07722-8B	Fixed 5" bariatric dual wheels	Y	1 ea
G07722G	Fixed 5" bariatric dual wheels	Y	3 pr/cs
G07723M	Footpiece, 5" wheels	Y	1 pr
MDS86615W	Footpiece, 3" wheels	Y	1 pr
G07820-8M	For all Guardian walkers	Y	1 pr
MDS86615SW5	Footpiece, 5" wheels	Y	2 pr/cs
G07732G	Glidebrakes	Y	3 pr/cs
MDS86615T	Rubber tips	Y	1 pr
G07900CS	Glide caps	Y	24/cs
G07904	Glide caps	Y	6/cs
MDS86615C	1" glide caps for walkers	Y	1 pr
MDS86615CW	1" glide caps for walkers	Y	6/cs
MDS86615SKI	Ski gliders for walkers	Y	6/cs
G07712	4-piece footpiece extension	Y	1 ea
MDS86615L	10" leg replacement	Y	1 pr
G07850M	Flip tray for Guardian walkers	Y	3/cs
G07850MB	Bariatric walker flip tray	Y	3/cs
G07702	Walker platform attachment for Guardian walkers	Y	1 ea
G07703M	Walker platform attachment youth	Y	1 ea
MDS86615P	Platform attachment for MDS aluminum walkers	Y	1 ea
G07741	Walker front carry pouch	Y	3/cs
G07715	Carrying basket	Y	2/cs
G07872	Walker seat	Y	2/cs
MDS86615K	Walker basket	Y	2/cs
MDS86615PB	Platform attachment for basic walkers	Y	1 ea
MDSWBALL	Walker ball glides	Y	6 pr/cs
MDS86615BAG	2-button walker bag	Y	4/cs
MDSRUBSKI	Walker glides	Y	6/cs
MDS86615W5T	5" long extension wheels	Y	1 pr



# Knee walkers

## Standard knee walkers

- Easy to maneuver with dual hand brakes for superb stopping control
- Can be used with right or left leg



MDS86000G4

Item No.	Description	Wt. capacity	Knee pad height	HCPCS	Retail Pkg.	Pkg.
MDS86000EB	8" wheels, Blue	300 lbs. (136.1 kg)	14.5"-21"	E0118	N	1 ea
MDS86000EJR	6" wheels, Red	350 lbs. (158.8 kg)	14"-17.5"	E0118	N	1 ea
MDS86000G4	8" wheels, Black	300 lbs. (136.1 kg)	17"-23"	E0118	N	1 ea
MDS86000SS	Seated scooter	300 lbs. (136.1 kg)	21"-28.5"	-	N	1 ea

## All-terrain knee walkers

Item No.	Description	Wt. capacity	Knee pad height	HCPCS	Retail Pkg.	Pkg.
ATKNEE	12" wheels, Black	350 lbs. (158.8 kg)	20"-26"	E0118	N	1 ea
ATKNEEJR	8" wheels, Black	350 lbs. (158.8 kg)	17.25"-24"	E0118	N	1 ea



ATKNEE

Folds down  
for easy  
portability



# Home care beds & accessories

## Basic beds

- Headspring supports recessed to resist pinching and injuries
- Foot assembly recessed to reduce chance of pinching
- 2-year warranty on electronics, limited lifetime warranty on the frame and welds



MDR107003E

Item No.	Description	Wt. capacity	Height range	HCPCS	Retail Pkg.	Pkg.
MDR107002E	Semi-electric basic bed	450 lbs. (204.1 kg)	15"-20"	E0295	N	1 ea
MDR107003E	Full-electric basic bed	450 lbs. (204.1 kg)	15"-20"	E0297	N	1 ea
MDR107003ELO	Full-electric basic low bed	450 lbs. (204.1 kg)	9.5"-22.5"	E0297	N	1 ea

## Medlite beds

- Premium lightweight home care beds
- 5-year warranty on electronics, limited lifetime warranty on the frame and welds (all Medlite beds)



MDR107003L

Item No.	Description	Wt. capacity	Height range	HCPCS	Retail Pkg.	Pkg.
MDR107002L	Medlite semi-electric bed	450 lbs. (204.1 kg)	14.5"-22.5"	E0295	N	1 ea
MDR107003L	Medlite full-electric bed	450 lbs. (204.1 kg)	14.5"-22.5"	E0297	N	1 ea
MDR107003LO	Medlite full-electric low bed	450 lbs. (204.1 kg)	9.5"-22.5"	E0297	N	1 ea

## Bariatric bed

- Powerful motors and sturdy steel frame for heavy-duty homecare bed needs
- 2-year warranty on electronics, limited lifetime warranty on the frame and welds



MDR107004

Item No.	Description	Wt. capacity	Height range	HCPCS	Retail Pkg.	Pkg.
MDR107004	Full-electric bariatric home care bed	600 lbs. (272.2 kg)	18"-26.5"	E0301	N	1 ea

### Side rails

- Half- and full-length Medline home care bed side rails
- 1-year warranty

Item No.	Description	HCPCS	Retail Pkg.	Pkg.
MDS89697	Clamp on half rail	E0305	N	1 pr
MDS89698N	Spring loaded half rail	E0305	N	1 pr
MDS89694N	Spring loaded full rail	E0310	N	1 pr
MDS89695N	Economy full rail	E0310	N	1 pr
MDS89697BAR	Bariatric half rail	E0305	N	1 pr



MDS89697



MDS89694

### Bed accessories

- Low bed conversion kits and transport aids for Medline home care beds

Item No.	Description	HCPCS	Retail Pkg.	Pkg.
MDSBEDCART	Bed cart kit for home care beds	-	N	1 ea
MDSLOWBASIC	Low bed conversion kit for MDR107003E	-	N	1 ea
MDSLOWBED	Low bed conversion kit for MDR107003L	-	N	1 ea



MDSLOWBED



MDSBEDCART

### Reacher

- Allows user to reach items both high and low

Item No.	Description	Handle material	HCPCS	Retail Pkg.	Pkg.
MDS86031RCE	31" reacher, Blue		-	N	6/cs



MDS86031RCE


### Bed assist bars

- Provides assistance getting in and out of bed
- Multiple hand positioning options to meet user's preference

Item No.	Description	Handle material	HCPCS	Retail Pkg.	Pkg.
MDS6800BA	Black	Foam	E0315	Y	3/cs
MDS6800BAM	M-shape handle, White	Foam	E0315	N	1 ea
MDS6800SWIV	SwivAssist bar, Black	Foam	E0315	N	3/cs
MDS6800SWIVRH	SwivAssist bar, Retail	Foam	E0315	Y	1 ea



MDS6800BA



MDS6800BAM



MDS6800SWIV



### Mattresses

- Cost-effective homecare bed mattresses available in three different materials
- Durable vinyl covers

Item No.	Description	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
MDR237827	Innerspring, 36" x 80" x 6"	250 lbs. (113.4 kg)	E0271	N	1 ea
MDR237827E	Fiber, 36" x 80" x 6"	250 lbs. (113.4 kg)	-	N	1 ea
MDR230981R	Foam, 36" x 80" x 6"	350 lbs. (158.8 kg)	E0272	N	1 ea
MDR23784	Innerspring 36" x 84" x 6"	250 lbs. (113.4 kg)	E0271	N	1 ea
MDR23784E	Fiber, 36" x 84" x 6"	250 lbs. (113.4 kg)	-	N	1 ea
MDR211892R	Foam, 42" x 84" x 6"	600 lbs. (272.2 kg)	E0272	N	1 ea

MDR237827



MDR237827E



MDR230981R



### Overbed tables

- Standard overbed tables with various features and material options
- Ideal for activities involving hemi, super-hemi, pediatric height wheelchairs and low beds

Item No.	Description	Wt. capacity	Height range	HCPCS	Retail Pkg.	Pkg.
MDS104015	Economy overbed table	50 lbs. (22.7 kg)	28"-45"	E0274	N	1 ea
MDS107015	Composite top overbed table	50 lbs. (22.7 kg)	30"-41"	E0274	N	1 ea
MDS108015	Fold-flat overbed table	50 lbs. (22.7 kg)	29.5"-41"	-	N	1 ea
MDS104950	Tilt-top overbed table	50 lbs. (22.7 kg)	30"-42"	E0274	N	1 ea

MDS104015



### Trapezes

- Assist with patient transfers, movement and positioning
- Free-standing and bed-mounted options

Item No.	Description	Wt. capacity	HCPCS	Retail Pkg.	Pkg.
MDS80615T	Standard trapeze bar	250 lbs. (113.4 kg)	E0910	N	1 ea
MDS80616T	Standard trapeze base	250 lbs. (113.4 kg)	A9900	N	1 ea
MDS80618T	Economy trapeze bar	250 lbs. (113.4 kg)	E0910	N	1 ea
MDS80619T	Economy trapeze base	250 lbs. (113.4 kg)	A9900	N	1 ea
MDS500TPZ	Lightweight bariatric trapeze	500 lbs. (226.8 kg)	E0912	N	1 ea
MDS1000TPZ	Bariatric trapeze	1000 lbs. (453.6 kg)	E0912	N	1 ea

MDS1000TPZ



# Canes

## Offset canes

- Single point canes
- With foam handles and wrist straps

Item No.	Description	Wt. capacity	Handle material	Height range	HCPCS	Retail Pkg.	Pkg.
MDS86420	Black	300 lbs. (136.1 kg)	Foam	29"-38"	E0100	Y	6/cs
MDS86420BRZW	Bronze	300 lbs. (136.1 kg)	Foam	29"-38"	E0100	Y	6/cs
MDS86420CHRW	Chrome	300 lbs. (136.1 kg)	Foam	29"-38"	E0100	Y	6/cs

## Quad canes

- Provides a secure grip with four feet
- Allows greater weight distribution for stability than single-point canes

Item No.	Description	Wt. capacity	Handle material	Height range	HCPCS	Retail Pkg.	Pkg.
MDS86222W	Small base, Black	300 lbs. (136.1 kg)	Foam	29"-38"	E0105	Y	2/cs
MDS86222CHR	Small base, Chrome	300 lbs. (136.1 kg)	Foam	29"-38"	E0105	Y	2/cs
MDS86228	Large base, Black	300 lbs. (136.1 kg)	Foam	29"-38"	E0105	Y	2/cs
MDS86228CHR	Large base, Chrome	300 lbs. (136.1 kg)	Foam	29"-38"	E0105	Y	2/cs

## Folding T-handle cane

- Folding capabilities make this cane perfect for daily travel
- Fit in a suitcase or carry-on luggage

Item No.	Description	Wt. capacity	Handle material	Height range	HCPCS	Retail Pkg.	Pkg.
MDS86430W	Black	250 lbs. (113.4 kg)	Vinyl	32"-36"	E0100	Y	6/cs



MDS86420BRZW



MDS86420



MDS86222W



MDS86430W



## Standard aluminum canes

- Standard handle lets users hook the cane over their arm while opening doors or performing other tasks
- Foam grip provides a comfortable secure grip

Item No.	Description	Wt. capacity	Handle material	Height range	HCPCS	Retail Pkg.	Pkg.
MDS86415	Black	300 lbs. (136.1 kg)	Foam	30"-38"	E0100	Y	6/cs
MDS86415CHR	Chrome	300 lbs. (136.1 kg)	Foam	30"-38"	E0100	Y	6/cs
G05370	Chrome	300 lbs. (136.1 kg)	Vinyl	30"-39"	-	Y	6/cs

## Bariatric canes

- All steel construction with dual-tubing insert for added strength and security

Item No.	Description	Wt. capacity	Handle material	Height range	HCPCS	Retail Pkg.	Pkg.
MDS86420XW	Offset handle, Black	500 lbs. (226.8 kg)	Foam	29"-38"	E0100	Y	1 ea
MDS86420XWCHR	Offset handle, Chrome	500 lbs. (226.8 kg)	Foam	29"-38"	E0100	Y	1 ea
MDS86222XW	Small base quad, Black	500 lbs. (226.8 kg)	Foam	29"-38"	E0105	Y	1 ea
MDS86228XW	Large base quad, Black	500 lbs. (226.8 kg)	Foam	29"-38"	E0105	Y	1 ea

## Hybrid canes

- Gel grip for user comfort
- Folding capabilities make this cane perfect for travel
- Four-point tip allows cane to be self-standing

Item No.	Description	Wt. capacity	Handle material	Height range	HCPCS	Retail Pkg.	Pkg.
4UCANE	Aluminum, Black	350 lbs. (158.8 kg)	Gel grip	32"-36"	E0100	N	1 ea
4UCANEB	Aluminum, Blue	350 lbs. (158.8 kg)	Gel grip	32"-36"	E0100	N	1 ea

## Cane accessories

Item No.	Description	Retail Pkg.	Pkg.
MDS86425W	Quad cane tips	Y	6/cs
MDS86426W4	Cane tips	Y	4/cs
MDS4PTCANETIP	4-point cane tip	Y	4/cs



MDS86415CHR G05370 MDS86420XW MDS86228XW 4UCANEB 4UCANE

# Our complete selection will help you every step of the way.

**Medline is your trusted DME partner**, with a wide range of premium products that enhance comfort and mobility, so patients can live their best life.



**We make  
healthcare  
run better™**

**Medline Industries, LP**  
Three Lakes Drive, Northfield, IL 60093  
Medline United States | 1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

**Medline Canada**  
1-800-268-2848 | medline.ca | canada@medline.com

**Medline México**  
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    