

WHEELCHAIR QUICK REFERENCE GUIDE

This comparison guide guide helps you navigate the various options and find the one that's right for you.

Product	Product Weight Capacity				Seat Width				Seat Depth			Arm Styles			Armpad Styles			Legrest/Footrest				Frame Finish				Upholstery Type		HCPCS Code
			16"	18"	20	, 22	" 24"	26"	28"	30"		Fixed	Removable	Removable Swing Back	Full-Length	Desk-Length	Desk-Length Ht. Adjustable	Fixed Footrest	Swing-Away Footrest	Elevating Legrest	Silver Vein Hammertone	Gray Powder Coated	Black Powder Coated	Blue Powder Coated	Chrome Plated	Nylon	Vinyl	
K1 Basic	38-lb.	300-lb.	•	•	•						16"	•	•		•	•			•	•	•					•		K0001 or K0002
K2 Basic	38-lb.	300-lb.	•	•	•						16"	•	•		•	•			•	•	•						•	K0001 or K0002
K3 Basic	35-lb.	300-lb.	•	•	•						16"		•			•			•	•	•					•		K0003
K3 Basic Plus	35-lb.	300-lb.	•	•	•						16"			•			•		•	•	•					•		K0003
K4 Basic*	33-lb.	300-lb.	•	•	•						16"/18"			•		•			•	•		•				•		K0004
K4 Lightweight*	31.5-lb.	300-lb.; 350-lb.	•	•	•						16"/18"			•		•			•	•			•			•		K0004
K4 Plus*	31.5-lb.	300-lb.; 350-lb.	•	•	•						16"			•			•		•	•			•			•		K0004
Excel 2000	35-lb.	300-lb.	•	•	•						16"	•	•		•	•		•	•	•					•		•	K0001 or K0002
Extra Wide	58-lb.	500-lb.			•		•				18"		•		•	•			•	•					•		•	K0007
Hybrid 2	33-lb.	300-lb.	•	•							16"		•			•			•	•				•		•		K0004
Recliner	38-lb.	300-lb.	•	•	•						17"		•			•				•					•		•	K0003 or K0007
Shuttle Bariatric	85-lb.	850-lb.					•	•	•	•	20"		•		•					•	•					•		K0007

All our wheelchairs come with the following warranty: limited lifetime on the frame, one year on the components.

*Larger seat depths are available.



K1 Basic Upholstery Color: Black



K2 Basic Upholstery Color: Black



K3 Basic Upholstery Color: Black



K3 Basic Plus Upholstery Color: Black



K4 Basic Upholstery Color: Black



K4 Lightweight Upholstery Color: Black





Excel 2000 Upholstery Color: Black/Navy



Extra Wide Upholstery Color: Navy



Hybrid 2 Upholstery Color: Black



Recliner Upholstery Color: Black



Shuttle Bariatric Upholstery Color: Black



K4 Plus Upholstery Color: Black

HCPCS DISCLAIMER: HCPCS codes and home health consolidated billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. The customer must address all coverage and reimbursement issues (including correctness and accuracy of codes) with their individual payers. It is the customer's responsibility to ensure the accuracy and appropriateness of each claim the customer submits, in accordance with all applicable payer requirements.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. ©2016 Medline Industries, Inc. Medline and Restore are registered trademarks of Medline Industries, Inc. MKT1667937 / e16134 / 35