

EDGE 3 *Stretto*[®]

EDGE 3 STRETTO[®]
shown with iLevel[®] power seat elevation,
Q-Logic 3 environmental controls, and
Stealth[®] TRU-Comfort 2 cushion.



The Edge 3 Stretto[®] features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12.5" or 14" drive wheels and optional 12" of iLevel[®] power adjustable seat height.



Narrowest
Base in
America

Q
QUANTUM[®]
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	14" or 12.5" pneumatic or flat free tires
CASTER WHEELS	Front: 5" flat free Rear: 6" flat free
MAXIMUM SPEED ¹	Up to 6.25 mph Up to 3.5 mph w/iLevel activated
GROUND CLEARANCE	2.9"
TURNING RADIUS ²	21"
OVERALL LENGTH ²	32.5" without front riggings
BASE WIDTH	21.75" with 14" drive wheels 20.47" with 12.5" drive wheels
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 12-20" Depth: 12-20"
TRU-Balance® 3 Power Tilt	Width: 12-20" Depth: 12-20"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-20" Depth: 12-20"
TRU-Balance® Flex Seating	Width: 15-20" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	16-18"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"
TRU-Balance® Flex Seating	16.5-18.25"
MANUAL RECLINE ⁶	Yes (TRU-Balance® Flex only)
BATTERY SIZE ³	40 AH (2) or 55 AH (2)
BATTERY WEIGHT ⁷	40 AH: 31 lbs. each 55 AH: 40 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel ⁶
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	126.4 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848) Single Power (K0856) Multiple Power (K0861)

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.

7) Battery weight may vary +/- 3 lbs. based on manufacturer.

8) All Quantum Rehab® units can be ordered with specialty controls.

9) Call for availability.

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- 19 bold, contemporary color choices
- Made in America of foreign and domestic parts


Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogy vent tray⁹
- Compact vent tray⁹
- 12" x 12" seat sizes with iLevel (SOF only)
- WC19 crash-test approved occupied transit kit
- Accu-Trac Advanced Tracking Technology

Colors



Hydrographics



QUANTUM®

#1 FOR REHAB POWER

(US) 866-800-2002

| (Canada) 888-570-1113

| quantumrehab.com

