



Q6 EDGE® Z
shown in *Back in Black*
with iLevel® technology

Q6 edge® Z

The Q6 Edge® Z utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6® design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge® Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability, and is available with our iLevel® technology.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Run flat Optional: 14" Pneumatic 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Standard: Up to 6 mph Optional: Up to 8 mph
GROUND CLEARANCE	3" at frame base
TURNING RADIUS ²	19.75"
OVERALL LENGTH ²	35.6"
BASE WIDTH	26.25"
SEATING SIZES	
SYNERGY ®	Width: 10-24" Depth: 10-24"
TRU-BALANCE ® 3 STATIC	Width: 12-24" Depth: 12-24"
TRU-BALANCE ® 3 POWER TILT	Width: 12-24" Depth: 14-24"
TRU-BALANCE ® 3 LIFT	Width: 12-24" Depth: 12-24"
SEAT-TO-FLOOR HEIGHTS	
SYNERGY ®/TRU-BALANCE® 3 STATIC ⁵	16.5-18.5"
TRU-BALANCE ® 3 POWER TILT ⁵	18.5-20.5"
TRU-BALANCE ® 3 iLEVEL® & TILT ⁵	18.5-20.5"
MANUAL RECLINE	No
BATTERY SIZE ³	Group 24
BATTERY WEIGHT ⁷	52 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE+ 75A Q-Logic NE
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole Mid-Wheel 6®
WEIGHT CAPACITY	300 lbs.
BASE WEIGHT	143.5 lbs.

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Part of Community Use Package with 6 mph motors. Increases width by 1.5".
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) Seat backs for units without TRU-Balance® 3 seating are only available in Back in Black.
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
 8) All Quantum Rehab® units can be ordered with specialty controls.
 9) Call for availability

Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 3 Power Positioning Systems
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- iLevel® option, featuring patented Extra Stability Technology™, enables up to 12" of seat elevation while at walking speed (up to 3.5 mph)
- Standard switch-operated LED fender lights

Options

- iLevel® power adjustable seat height
- TRU-Balance® 3 Seating
- Q-Logic 3 Drive Controls
- 8 mph motor package
- Accu-Trac Advanced Tracking Technology⁹
- Trilogy vent tray⁹
- Compact vent tray⁹
- 4" wide flat-free or pneumatic knobby extended traction tires^{9,4}

Colors

Available on Q-Logic 3 controllers and units with TRU-Balance® 3 (base and back)⁶

